



WARRANTY

TECNIPLAST H-TEMP Polysulfone (PSU) caging is unconditionally guaranteed for 3 years from the stamped date.

During the warranty period Tecniplast liability shall be limited to the replacement of all severely warped, discoloured or cracked products not caused by improper use.

All defective products must be returned to Tecniplast spa, Via I Maggio n.6 - 21020 Buguggiate - VA (Italy).



Tecniplast S.p.a.

Via 1° Maggio, 6 • 21020 Buguggiate - VA • Italy
Tel. +39 0332 80 97 11 • Fax +39 0332 45 83 15
www.tecniplast.it • E-mail: tecnicom@tecniplast.it



Hereby follows a non exhaustive list of improper usage examples:

- Washing with highly alkaline Ph detergents (particularly prolonged exposure to strong alkali in a hot environment).
- Washing with any ester-based chemicals.
- Rinse aids.
- Autoclaving to temperatures higher than 134°C/273°F.
- Autoclaving in non-suitable plastic bags (sealed, no high steam permeability, wrong dimensions).

Please contact our local distributor in case you need additional information.



Tecniplast S.p.a.

Via 1° Maggio, 6 • 21020 Buguggiate - VA • Italy
Tel. +39 0332 80 97 11 • Fax +39 0332 45 83 15
www.tecniplast.it • E-mail: tecnicom@tecniplast.it